


US-FDA APPROVED REGEN LAB BCT PRP KIT





US-FDA APPROVED
regenlab
PRP




Platelet Rich Plasma,

- **Autologous Platelet Rich Plasma (A-PRP) is plasma with a concentration of platelets above the blood baseline.**
- **A-PRP is developed from autologous blood, it is therefore inherently safe and free of transmissible diseases such as HIV and hepatitis.**
- **Within A-PRP, the concentration of platelets delivers an increased number of growth factors in their biologically determined ratios. This distinguishes PRP from recombinant (not native) growth factors.**
- **The plasma fibrinogen will serve to form the natural scaffold for tissue reconstruction.**

Growth Factors in Platelets:


Factor	Name	Principal Source	Effects
PDGF aa PDGF bb PDGF ab	Platelet derived growth factors	Activated thrombocytes	Mitogenes of mesenchymal stem cells promote the synthesis of the extracellular matrix
TGF-alpha TGF-beta	Transforming Growth Factors	Activated thrombocytes	Stimulation of DNA synthesis, proliferation of various types of cells. Favorises the synthesis of collagen
IGF-I IGF-II	Insulin-like Growth Factors	Activated thrombocytes	Stimulates the proliferation and differentiation of osteoblasts
EGF	Epidermal Growth Factor	Activated thrombocytes	Stimulates proliferation and differentiation of epidermis cells, co-stimulating angiogenesis
VEGF	Vascular Endothelial Growth Factor	Leucosytes & Endothelial cells	Stimulate angiogenesis & chemo-attraction of osteoblasts

REGENERATIVE MEDICINE: REPAIRING & RESTRUCTURING:



PATIENT AND DOCTOR EXPECTATIONS FOR A NEW TREATMENT, IMPROVEMENT IN:

- **Tolerance**
- **Skin Revitalisation**
- **Product Durability**
- **Safety**



**REGEN LAB PRP
AUTOLOGOUS PLATELET-RICH PLASMA**

REGEN BCT Platelet- Rich Plasma:

- **Autologous - safe, biocompatible**
- **Increases tissue vascularity through angiogenesis**
- **Chemotactic for stem cell, macrophages, and fibroblasts**
- **Increases the rate of epithelial, collagen and granulation tissue production**
- **Antimicrobial effect**
- **Surgical hemostatic agent for wound hemostasis and lymphatic sealing.**

For preparation of PRP

WHY REGEN LAB PRP KIT

ONLY?

BECAUSE:-

- **It is the only kit which is approved by US-FDA for collection PRP**
- **SIMPLICITY** in methodology
- **LOW BLOOD VOLUME REQUIRED**
- **OPTIMUM PLATELET RECOVERY**
- **CONSISTENT ISOLATION OF A-PRP**
- **RED BLOOD CELL DEPLETION**
- **Regen Lab has more then 14 years clinical experience in aesthetics clinicians**
- **PRP has been proven for various clinical trials.**

- Regen Kit® is designed to prepare Autologous Platelet Rich Plasma by centrifugation of peripheral blood in a vacuum tube, using a thixotropic gel for cell separation and citrate as anticoagulant.
- The system is simple to use, works in a closed system and produces 4 to 5 ml of Platelet Rich Plasma from 8 ml of whole blood in less than 10 minutes.
- The gel in the Regen BCT tubes specifically enables the recovery of more than 95% of the platelets and isolates them from red blood cells.

The table below shows the results of a study on A-PRP obtained from 60 individuals (34 men and 26 women, age 18 to 76, mean age 35±14), processed according to the Regen Kit® protocol.

PRP volume	Platelet Recovery	Red blood Cell remnant	Platelet Concentration Factor
4 to 5 ml	> 95 %	< 2 %	1.7 to 2.5 x

ADVANTAGES

SIMPLICITY

Only few straightforward steps are required to prepare 4 to 5 ml of Platelet-Rich Plasma.

LOW BLOOD VOLUME REQUIRED

RegenKit® is designed to prepare Autologous Platelet Rich Plasma from only 8 ml of blood.

OPTIMUM PLATELET RECOVERY

With Regen Kit® you recover almost 100% of the platelets from the whole blood. Platelets are concentrated in 1 to 5 ml of autologous plasma (depending on the required platelet concentration).

CONSISTENT ISOLATION OF A-PRP

Regen Kit® uses a physical barrier (thixotropic gel) to separate the red blood cells from the plasma and the cellular sediment. The preparation is safe and consistent.

RED BLOOD CELL DEPLETION

The level of Red Blood Cells is drastically reduced, with more than 98% of RBCs removed.

- **Simple and fast process. PRP is ready with only few steps in less than 10 minutes.**
- **Only 8ml of whole blood are required to prepare 4 To 5 ml of PRP**
- **Double sterile blister. Regen PRP is the only PRP approved to be used in an operating room**
- **The Regen technology effectively isolates the platelet concentrate from the RBCs**

- **PRP has been proven for various clinical trials,**
- **Regen Lab has more than 12 years clinical experience in aesthetics clinicians.**
- **Regen offers proven detailed clinical protocols, clinical studies, a dedicated YouTube channel with educational videos, as well as in house clinical specialists and consultant doctors**

Highest Safety Measures:

- **Pyrogen free.**
- **Sterility ensures the absence of viable living bacteria. However, the act of sterilizing a contaminated vial can actually result in the release and deposit of pyrogens, as sterilization destroys bacteria, leading to bacterial cell lysis and release endotoxins and exotoxins.**


International Approvals:

- **CE Mark**
- **US-FDA**
- **Health Canada**
- **K FDA**
- **Australia, Japan, Kuwait, Saudi Arabia, Iran, UAE, Singapore, Malaysia, Indonesia, Hong Kong, South Africa etc**


A-PRP® PREPARATION PROTOCOL




Centrifuge @ 3400 rpm
(1500g) for 5 minutes
in Regen Centrifuge or
6 minutes in Remi-R8 C
centrifuge



HOMOGENIZE




Regen BCT TUBE Before & After centrifugation:




**Most efficient PRP
system, loosing only 5%
of platelets during
processing.....**


Cell survival:



Déplétion of 99.7% of the hématies:



Platelet réponse to induced activation:



Platelets: Are more really better?

- PRP with near physiological concentration show better results for wound healing than too highly concentrated PRP:
- “Those methods with lower concentrations of platelets – 1 to 3 times baseline – showed more robust healing rates than those with higher concentrations of 3 to 8 times baseline.”
- Rappl LM et al. Effect of platelet-rich plasma gel in a physiologically relevant platelet concentration on wounds in persons with spinal cord injury. Int Wound J 2011; 8:18.7–195.

- **Animal studies show negative effect of too high platelet concentration:**
- **“The use of highly concentrated platelet preparations appeared to have an inhibitory influence on osteoblasts activity. Possible reasons could be unwanted inhibitory and cytotoxic effects of growth factors at such high concentrations.”**
- **Weibrich G. et al., Effect of platelet concentration in platelet-rich plasma on peri-implant bone regeneration. Bone 2004; 34:665-671.**

- “PRP might exert positive effects on intestinal anastomotic healing in a dose-dependent manner up to a certain level, but adverse effects occur when it is highly concentrated.”

Yamaguchi R. al., Effects of Platelet-Rich Plasma on Intestinal Anastomotic Healing in Rats: PRP Concentration is a Key Factor. J Surg Res. 2012 Apr;173(2):258-66.

- In vitro studies reach the same conclusion:

“Optimal results were observed at a platelet concentration of 1.7 X”,
“Increased concentrations resulted in a reduction in proliferation and a suboptimal effect on osteoblasts function.”

Graziani F. et al., The in vitro effect of different PRP concentrations on osteoblasts and fibroblasts. Clin. Oral. Impl. Res. 17, 2006; 212–219.

Absolute contraindications are:

- **Marked thrombocytopenia**
- **Blood diseases**
- **Hypofibrinogenemia**
- **Hemodynamic instability**
- **Anticoagulant or antiaggregant treatment**
- **Oncological diseases (specially, skeletal system and blood)**
- **Sepsis**
- **Acute and chronic infectious diseases**
- **Autoimmune diseases, for example, lupus erythematosus, etc.**

Relative contraindications:

- **Chronic liver diseases in the exacerbation phase**
- **Chronic intoxication against the background of long-term use of alcohol, addictive drugs, an administration of potent medications**
- **Use of no steroidal anti-inflammatory drugs in little less than 2 days before drawing the blood, an injection of cortico-steroids in little less than 2 weeks before the procedure**
- **Pregnancy**
- **Inflammatory skin diseases, chronic dermatosis in the exacerbation phase, for example, psoriasis, atopic eczema, etc.**

Literature on PRP:

**5000 search results on Pubmed
experiments 1957-1998**

Lab & animal

Dental publications 1998 Marx,

Wound healing 2001 Powell

Cardiac publications 2005 Englert

Orthopedic publication 2006 Anitua

Aesthetic publications 2007 Gondinet

- **Carter, M. J., C. P. Fylling, et al. (2011). "Use of platelet rich plasma gel on wound healing: a systematic review and meta-analysis." Eplasty 11: e38.**
- **Cervelli, V., P. Gentile, et al. (2009). "Application of platelet-rich plasma in plastic surgery: clinical and in vitro evaluation." Tissue Eng Part C Methods 15(4): 625-634.**
- **Cervelli, V., F. Nicoli, et al. (2012). "Treatment of traumatic scars using fat grafts mixed with platelet-rich plasma, and resurfacing of skin with the 1540 nm nonablative laser." Clin Exp Dermatol 37(1): 55-61.**
- **Cho, J. M., Y. H. Lee, et al. (2010). "Effect of platelet-rich plasma on ultraviolet b-induced skin wrinkles in nude mice." J Plast Reconstr Aesthet Surg.**
- **Farrior, E. and K. Ladner (2011). "Platelet gels and hemostasis in facial plastic surgery." Facial Plast Surg 27(4): 308-314.**
- **Marx, R. E. (2004). "Platelet-rich plasma: evidence to support its use." J Oral Maxillofac Surg 62(4): 489-496.**
- **Mazzucco, L., V. Balbo, et al. (2008). "Platelet-rich plasma and platelet gel preparation using Plateltex." Vox Sang 94(3): 202-208.**
- **Redaelli (2010). "Face and Neck Revitalization With Platelet-rich Plasma." Journal of Drugs in Dermatology 9(5): 466-472.**



Before and after pictures SKIN REJUVENATION



Nicole Belle RN



Nicole Belle RN



Nicole Belle RN




Before and after pictures





Before and after pictures



Before and after pictures



Before and after pictures



Before and after pictures


Photos courtesy of Nicole Bell



Before treatment



After treatment results
Regenerative Medicine Lower Face and Neck
Improved Texture, Tone & Tighten




Before and after pictures
AGA

Hair Regeneration

Protocol
Standard BCT x 2
Nutritional Supplementation
















The Australian Experience with Regen PRP

Neck Rejuvenation in 6 weeks



**A 50 year old woman's birthday present to herself
- It exceeded her expectations!**



**Protocol
2 BCT tubes every 2 weeks, with RF skin tightening treatment in between
Activation in vivo**

Neck Rejuvenation in 6 weeks

**Results from superficial micro blebs (1 bevel deep) through
necklace lines
Result after 6 weeks**



Peri-orbital Rejuvenation

The thinning skin of the ageing eyelids are difficult to treat with traditional modalities such as dermal fillers or fractional laser

Stand-alone PRP represents an effective and safe method of rejuvenation of the eyes in patients from their late 20s and upwards

Protocol – Session x 1; Activation in vivo; Micro bleb



Peri-orbital Rejuvenation



Stand-alone PRP has resulted in improved skin texture, elasticity, micro-vascularisation and skin thickening after just a single treatment for this patient

Protocol
BCT x 1
Sessions x 1
Activation with calcium
Micro bleb



Peri-orbital Rejuvenation



Fresher in 6 weeks

**Alternating weekly
treatments of Regen and RF
skin tightening**

Protocol


**No Botulinum Toxin used in
this period**

BCT x 1 around eyes

**Tangential micro blebs with
30G short needle**

Activation with Calcium





Facial Rejuvenation Combination Treatment



**Combination treatment
PRP with Portrait Skin Resurfacing.**
This patient demonstrates overall improvement in skin texture, fine lines and pigmentation, with accelerated healing



Protocol
2 x BCT tubes whole face
Activation with calcium
Portrait Resurfacing

Facial Rejuvenation



Stand-alone PRP

Note the improvement in fine lines, pore size and texture



Protocol
3 sessions
Activation in vivo

Acne Scar Improvement

Regen injected superficially and undermining scars in combination with microneedling

Protocol: BCT x 2, 3 sessions, Activation in vivo



Neck Rejuvenation



Neck rejuvenation has long been the province of the surgeon.

PRP in combination with RF and HA dermal fillers can improve the skin of the neck, often with greater effectiveness than seen on facial skin



Protocol

BCT x 4 per treatment

Activation calcium

RF skin tightening, HA fine lines

Period of treatment: June 2009 – June 2013


Body (décolletage) Rejuvenation

The décolletage is a neglected area in Cosmetic Medicine. Depending on environmental (ultraviolet) exposure, signs of ageing will begin appearing in the 20s. The focus is on prevention but the décolletage responds very quickly to PRP and less treatments are generally required
BCT x 2, 1 session, long needle



Body (limb) Rejuvenation

The ageing skin of the limbs has again been a neglected area. PRP in combination with RF skin tightening can significantly improve the appearance of the ageing arm (“Tuck-shop lady arms”)



- **You Tube video link for preparation of PRP:-**

https://www.youtube.com/watch?v=S3sJkP3q1h4&list=PLusMjgkfTtg3zXvaCvQLesZP_s2ECro04&index=2

FOR MORE DETAILS CONTACT:-



MR.PAWAN WAGRULKAR
PRODUCT HEAD,
GEOSMATIC C&C PVT. LTD.

EMAIL ID:- pawan@geosmatic.com /
knowledge@geosmatic.com

Mob. No.:- 077740 77252

Thank you.